

Rehabilitation Centres

Illicit substance and alcohol abuse in the workplace is becoming an increasingly common problem nationwide. Being under the influence of drugs or alcohol in a workplace environment can lead to inconsistent work quality, poor concentration and carelessness. It can be extremely dangerous, especially in a company where employees drive, work directly with consumers, and operate machinery. With far-reaching implications, up to 20 percent of work-related fatalities test positive for drugs or alcohol. The drugs at the root of this issue include cocaine, heroin, methamphetamine and prescription drugs.

Where do we fit in?

Each day millions of Americans walk around with an addiction and only a small percentage of individuals pursue treatment. In rehab drug testing is used as a tool rather than a treatment method. Drug testing eliminates the ability of addicts to manipulate the detox process, keeping individuals accountable, honest and drug-free during treatment. It is important to have a clear understanding of the drugs in the patient's system before facilitating detox or administering adjunct medications during treatment.

Benefits



Simple Process

1. Add sample to cartridge
2. Load sample and tip cartridges
3. Press start to get results



Multiple Matrices

Applicable for urine, blood and oral fluid enables less invasive and quick sample collection



Rapid Results

Qualitative results are available in under 20 minutes from one single employee sample.



Small sample volume

Only a small sample volume required, leaving more for confirmatory testing.

Drugs Detected on the Evidence MultiSTAT

| Drug Groups | | |
|---------------------------------------|-------------------------|---------------------------------|
| α -PVP (Flakka) | Cannabinoids (THC) | Opiates |
| AB-CHMINACA | Creatinine | Oxycodone |
| AB-PINACA | Ethyl Glucuronide (EtG) | Phencyclidine (PCP) |
| Amphetamine | Fentanyl | Propoxyphene (PPX) |
| Barbiturates | JWH-018 | Pregabalin |
| Benzodiazepines I (Oxazepam) | Ketamine | Tramadol |
| Benzodiazepines II (Lorazepam) | LSD | Tricyclic Antidepressants (TCA) |
| Benzoyllecgonine (Cocaine Metabolite) | Methadone | URI44 / XLR11 |
| Buprenorphine | Methamphetamine | 6-MAM |

*Our tests are validated for Intercept i2 and NeoSal™ collection devices.